

MEDICO CENTRIFUGE

Genesys offers a variety model of medico centrifuge which is ideal for routine centrifuging test of medical practitioners, hospitals, colleges, institutions, and industrial laboratories etc. Fitted with resiliently mounted motor for vibration free performance. A specially wound 5 step transformer is used for speed regulation. The highly sophisticated elegant body of the equipment is made with powder coated mild steel which gives excellent protection.



TECHNICAL SPECIFICATIONS

- ▮ Centrifuge Features : Swing out head, outer body made of Mild Steel powder coated.
- ▮ Maximum speed : 4000 RPM
- ▮ Maximum RCF : 1880 x g